


Mephenytoin-D3

Chemical Properties

CAS No. :	1185101-86-3
Formula:	C ₁₂ H ₁₄ N ₂ O ₂
Molecular Weight:	221.27
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Mephenytoin-D3 is the deuterated form of Mephenytoin (T21425) and can be used for isotope tracing.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.5194 mL	22.5968 mL	45.1937 mL
5 mM	0.9039 mL	4.5194 mL	9.0387 mL
10 mM	0.4519 mL	2.2597 mL	4.5194 mL
50 mM	0.0904 mL	0.4519 mL	0.9039 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481