

Glycidyl methacrylate

Chemical Properties


CAS No. : 106-91-2

Formula: C₇H₁₀O₃

Molecular Weight: 142.15

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Glycidyl methacrylate is an ester of methacrylic acid and glycidol. Glycidyl methacrylate is a common monomer used in the creation of epoxy resins.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.0348 mL	35.1741 mL	70.3482 mL
5 mM	1.407 mL	7.0348 mL	14.0696 mL
10 mM	0.7035 mL	3.5174 mL	7.0348 mL
50 mM	0.1407 mL	0.7035 mL	1.407 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

- Ribeiro LF, Masini JC. Complexing porous polymer monoliths for online solid-phase extraction of metals in sequential injection analysis with electrochemical detection. *Talanta*. 2018 Aug 1;185:387-395. doi: 10.1016/j.talanta.2018.03.099. Epub 2018 Mar 31. PubMed PMID: 29759217.
- Obbed MS, Aqel A, Al Othman ZA, Badjah-Hadj-Ahmed AY. Preparation, characterization and application of polymethacrylate-based monolithic columns for fast and efficient separation of alkanes, alcohols, alkylbenzenes and isomeric mixtures by gas chromatography. *J Chromatogr A*. 2018 Jun 22;1555:89-99. doi: 10.1016/j.chroma.2018.04.045. Epub 2018 Apr 23. PubMed PMID: 29724647.
- Wu R, Ding X, Qi Y, Zeng Q, Wu YW, Yu B, Xu FJ. Flexible Cationic Nanoparticles with Photosensitizer Cores for Multifunctional Biomedical Applications. *Small*. 2018 May 2:e1800201. doi: 10.1002/sml.201800201. [Epub ahead of print] PubMed PMID: 29717807.
- Kim SH, Yeon YK, Lee JM, Chao JR, Lee YJ, Seo YB, Sultan MT, Lee OJ, Lee JS, Yoon SI, Hong IS, Khang G, Lee SJ, Yoo JJ, Park CH. Precisely printable and biocompatible silk fibroin bioink for digital light processing 3D printing. *Nat Commun*. 2018 Apr 24;9(1):1620. doi: 10.1038/s41467-018-03759-y. PubMed PMID: 29693652; PubMed Central PMCID: PMC5915392.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481