


2H,6H-Benzo[1,2-b:5,4-b']dipyran-2-one, 7,8-dihydro-8,8-dimethyl-

## Chemical Properties

CAS No. : 21074-63-5  
 Formula: C<sub>14</sub>H<sub>14</sub>O<sub>3</sub>  
 Molecular Weight: 230.26  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	2H,6H-Benzo[1,2-b:5,4-b']dipyran-2-one, 7,8-dihydro-8,8-dimethyl- is a high-purity biochemical reagent, suitable for use as a biological material or organic synthesis intermediate, and widely used in research and experiments in the life sciences.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3429 mL	21.7146 mL	43.4292 mL
5 mM	0.8686 mL	4.3429 mL	8.6858 mL
10 mM	0.4343 mL	2.1715 mL	4.3429 mL
50 mM	0.0869 mL	0.4343 mL	0.8686 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481