


## Propanoic-13C3 sodium

## Chemical Properties

CAS No. :	152571-51-2
Formula:	C <sub>3</sub> H <sub>5</sub> NaO <sub>2</sub>
Molecular Weight:	99.04
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Propanoic-13C3 sodium is a 13C-labeled Sodium propionate (TMSM-2098).
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	10.0969 mL	50.4847 mL	100.9693 mL
5 mM	2.0194 mL	10.0969 mL	20.1939 mL
10 mM	1.0097 mL	5.0485 mL	10.0969 mL
50 mM	0.2019 mL	1.0097 mL	2.0194 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481