

PEG6000

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PEG6000 is used as a pharmaceutical excipient, such as in ointment bases and lubricants. Pharmaceutical excipients, or adjuvants, are chemical substances used in the pharmaceutical process excluding active drug components. They generally refer to the inactive components in drug formulations that enhance the stability, solubility, and manufacturability of the drug. Excipients also impact the absorption, distribution, metabolism, and elimination (ADME) processes of combined medications.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481