

OR51E1 Human Pre-designed siRNA Set A

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | "OR51E1 Human Pre-designed siRNA Set A consists of three pairs of siRNA targeting different regions of the OR51E1 (Human) gene. It also includes a negative control, a FAM-labeled negative control, and a positive control." |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481