


## L-Methionine-D3

## Chemical Properties

CAS No. :	13010-53-2
Formula:	C5H11NO2S
Molecular Weight:	152.23
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	L-Methionine-D3 is the deuterated form of L-Methionine (TMIM-00050). L-Methionine, the L-isomer of Methionine, is an essential amino acid. Known for its strong liver detoxifying properties, Methionine acts as a hepatic protectant.
Targets(IC50)	Endogenous Metabolite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.569 mL	32.845 mL	65.6901 mL
5 mM	1.3138 mL	6.569 mL	13.138 mL
10 mM	0.6569 mL	3.2845 mL	6.569 mL
50 mM	0.1314 mL	0.6569 mL	1.3138 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481