

Chitosan ($\geq 90\%$ deacetylated, Low viscosity, $< 200\text{mPa}\cdot\text{s}$)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Chitosan (Deacetylated chitin) is a biochemical reagent that serves as a biomaterial or organic compound, useful in life science research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481