

## Vericiguat-D6

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Vericiguat-D6 is the deuterium-labeled version of Vericiguat. Vericiguat (T17226) (BAY1021189) functions as a potent and orally active stimulator of soluble guanylate cyclase.
Targets(IC50)	Guanylate cyclase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481