

Boc-Asp(OBzl)-Pro-Arg-AMC acetate

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from moisture**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Boc-Asp(OBzl)-Pro-Arg-AMC acetate is a fluorescent substrate used for assessing protease activity. Upon hydrolysis, it liberates the fluorescent product 7-amino-4-methylcoumarin (AMC), which emits a fluorescent signal under ultraviolet light.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481