

Atropine-D5

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Atropine-D5 is the deuterated form of Atropine sulfate monohydrate (T0375L) . Atropine (Tropine tropate) sulfate monohydrate serves as a broad-spectrum, competitive antagonist of muscarinic acetylcholine receptors (mAChR) and exhibits anti-myopia properties.
Targets(IC50)	AChR, Autophagy

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481