

Alpha-Estradiol-D3

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Alpha-Estradiol-D3 is the deuterated form of Alpha-Estradiol. Alpha-Estradiol (T3378) is a female hormone and a 5 α -reductase inhibitor, with potential applications in research addressing androgenic alopecia.
Targets(IC50)	Reductase,Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481