

DSPE-PEG-DBCO ammonium

Chemical Properties

CAS No. : 2052955-83-4

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DSPE-PEG-DBCO ammonium is the ammonium salt form of DSPE-PEG-DBCO. It is utilized in copper-free click chemistry by coupling with azide-functionalized peptide ligands through SPAAC conjugation reactions. This compound is employed in drug delivery systems and nanoparticle research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481