

Cilnidipine-D7

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cilnidipine-D7 is the deuterated form of Cilnidipine. As a long-acting, second-generation dihydropyridine Ca ²⁺ channel blocker, Cilnidipine (T0388) effectively targets both L-type and N-type Ca ²⁺ channels and exhibits antihypertensive properties.
Targets(IC50)	Calcium Channel

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481