

## Tofacitinib-13C3

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

|               |  |
|---------------|--|
| Description   | Tofacitinib-13C3 is a carbon-13 labeled variant of Tofacitinib. Tofacitinib (T6321) acts as a JAK3/2/1 inhibitor, exhibiting IC50 values of 1 nM, 20 nM, and 112 nM, respectively. |
| Targets(IC50) | Apoptosis,JAK  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481