

LysoPC(14:0/0:0)-D7

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	LysoPC(14:0/0:0)-D7 is a deuterium-labeled version of LysoPC (14:0/0:0) (TQ0275). This compound is a lysophospholipid (LyP), characterized as a monoglycerophospholipid with a phosphocholine moiety occupying the glycerol substitution site. LysoPC (14:0/0:0) exhibits potent antispasmodic properties.
Targets(IC50)	Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481