

L-Kynurenine-13C4,15N-1

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	L-Kynurenine-13C4,15N-1 is L-Kynurenine (T4928) labeled with 13C and 15N. It is a metabolite of L-tryptophan and acts as an agonist of the aryl hydrocarbon receptor (aryl hydrocarbon receptor).
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481