

Ubrogepant-D5

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Ubrogepant-D5 is the deuterium-labeled version of Ubrogepant. It is an orally active and selective antagonist of the calcitonin gene-related peptide receptor (CGRP). Ubrogepant (T17194) exhibits high affinity for CGRP receptors in humans and rhesus monkeys, effectively blocking the cAMP response induced by α -CGRP stimulation. This compound is utilized in studies focused on acute migraine.
Targets(IC50)	CGRP Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481