

PEG600

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	PEG600 can serve as an ointment base and lubricant, among other pharmaceutical excipients. These excipients or adjuvants are chemical substances used in the drug manufacturing process that are not active medicinal ingredients. They typically enhance the stability, solubility, and manufacturability of pharmaceutical formulations, while also influencing the absorption, distribution, metabolism, and excretion (ADME) of combined medications.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481