


## Deoxycholic acid-D5

## Chemical Properties

CAS No. :	52840-14-9
Formula:	C <sub>24</sub> H <sub>40</sub> O <sub>4</sub>
Molecular Weight:	397.60
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Deoxycholic acid-D5 is the deuterated form of Deoxycholic acid (TMSM-2674). It activates the G protein-coupled bile acid receptor TGR5.
Targets(IC50)	Endogenous Metabolite, GPCR19

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5151 mL	12.5755 mL	25.1509 mL
5 mM	0.503 mL	2.5151 mL	5.0302 mL
10 mM	0.2515 mL	1.2575 mL	2.5151 mL
50 mM	0.0503 mL	0.2515 mL	0.503 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481