


Calcium phosphorylcholine-D9 chloride

Chemical Properties

CAS No. :	344299-43-0
Formula:	C ₅ H ₁₃ CaClNO ₄ P
Molecular Weight:	266.72
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Calcium phosphorylcholine-D9 chloride is the deuterated form of calcium phosphorylcholine chloride (T65394).
Targets(IC50)	Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7493 mL	18.7463 mL	37.4925 mL
5 mM	0.7499 mL	3.7493 mL	7.4985 mL
10 mM	0.3749 mL	1.8746 mL	3.7493 mL
50 mM	0.075 mL	0.3749 mL	0.7499 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481