


## L-Tyrosine-13C

## Chemical Properties

CAS No. :	110622-46-3
Formula:	C <sub>9</sub> H <sub>11</sub> NO <sub>3</sub>
Molecular Weight:	182.18
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	L-Tyrosine-13C is the carbon-13 labeled form of L-Tyrosine. L-Tyrosine (T0493) is a non-essential amino acid known for inhibiting the activity of citrate synthase in the posterior cortex.
Targets(IC50)	Endogenous Metabolite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4891 mL	27.4454 mL	54.8908 mL
5 mM	1.0978 mL	5.4891 mL	10.9782 mL
10 mM	0.5489 mL	2.7445 mL	5.4891 mL
50 mM	0.1098 mL	0.5489 mL	1.0978 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481