

## Pulegone-D8

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Pulegone-D8 is a deuterium-labeled form of Pulegone (TMSM-1801). Tetracycline is an orally active broad-spectrum antibiotic effective against various bacteria, including gram-positive, gram-negative, chlamydiae, mycoplasmas, and rickettsiae. It is used for infection research.
Targets(IC50)	Endogenous Metabolite,TRP/TRPV Channel

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481