

## SLC39A10 Human Pre-designed siRNA Set A

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

|               |   |
|---------------|---|
| Description   | SLC39A10 Human Pre-designed siRNA Set A consists of three pairs of siRNAs targeting various regions of the SLC39A10 (Human) gene. The set also includes a negative control, a FAM-labeled negative control, and a positive control. |
| Targets(IC50) | Others  |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481