

OR51E2 Human Pre-designed siRNA Set A

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	OR51E2 Human Pre-designed siRNA Set A consists of three pairs of siRNA that are specifically designed to target various regions of the OR51E2 (Human) gene. Additionally, the set includes a negative control, a FAM-labeled negative control, and a positive control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481