

## OR51E2 Human Pre-designed siRNA Set A

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	OR51E2 Human Pre-designed siRNA Set A consists of three pairs of siRNA that are specifically designed to target various regions of the OR51E2 (Human) gene. Additionally, the set includes a negative control, a FAM-labeled negative control, and a positive control.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481