

(R)-MPH-220

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	(R)-MPH-220 represents the R isomer of MPH-220, a selective orally active inhibitor of skeletal muscle myosin-2. MPH-220 facilitates muscle relaxation and serves as an antispasmodic agent, making it useful for research involving spasms and muscle stiffness.
Targets(IC50)	Myosin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481