

Glb1 Mouse Pre-designed siRNA Set A

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Glb1 Mouse Pre-designed siRNA Set A includes siRNAs targeting the Glb1 gene (Mouse), along with a negative control, a positive control, and a FAM-labeled negative control.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481