

## Isobutyryl-L-carnitine-D3 chloride

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Isobutyryl-L-carnitine-D3 chloride is the deuterated form of Isobutyryl-L-carnitine chloride (T73766), which is a product of acyl-CoA dehydrogenase.
Targets(IC50)	Endogenous Metabolite

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481