

Deoxycytidine triphosphate-15N3 dilithium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Deoxycytidine triphosphate-15N3 dilithium is a variant of deoxycytidine triphosphate labeled with 15N. dCTP is a nucleotide triphosphate utilized in DNA synthesis and is applicable in research areas such as real-time PCR, cDNA synthesis, and DNA sequencing.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Endogenous Metabolite,DNA/RNA Synthesis

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481