

## Teprasiran sodium

## Chemical Properties

CAS No. : 2169310-72-7

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Teprasiran sodium is a small interfering RNA (siRNA) that inhibits p53-mediated cell death in acute kidney injury.
Targets(IC50)	MDM-2/p53

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481