

Glycochenodeoxycholic acid-D7 sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Glycochenodeoxycholic acid-D7 sodium is the deuterated form of Glycochenodeoxycholic acid (sodium salt). Glycochenodeoxycholic acid sodium (T5259) salt, also known as Chenodeoxycholyglycine sodium salt, is a bile salt produced in the liver from chenodeoxycholic acid and glycine. It functions as a detergent to dissolve fats for absorption and is absorbed itself.
Targets(IC50)	Apoptosis,Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481