

CHRNA7 Human Pre-designed siRNA Set A

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | CHRNA7 Human Pre-designed siRNA Set A contains three pairs of siRNAs targeting various regions of the CHRNA7 (Human) gene, along with a negative control, a FAM-labeled negative control, and a positive control. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481