


2-Arachidonoylglycerol-D8

Chemical Properties

CAS No. :	924894-97-3
Formula:	C ₂₃ H ₃₈ O ₄
Molecular Weight:	386.60
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2-Arachidonoylglycerol-D8 is the deuterated form of 2-Arachidonoylglycerol. 2-Arachidonoylglycerol serves as an endogenous cannabinoid ligand within the central nervous system.
Targets(IC50)	Cannabinoid Receptor, Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5867 mL	12.9333 mL	25.8665 mL
5 mM	0.5173 mL	2.5867 mL	5.1733 mL
10 mM	0.2587 mL	1.2933 mL	2.5867 mL
50 mM	0.0517 mL	0.2587 mL	0.5173 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481