

PEG1500

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PEG1500 serves as an ointment base and lubricant, functioning as a pharmaceutical excipient. Pharmaceutical excipients or adjuvants are defined as chemical substances used in the drug formulation process apart from the active pharmaceutical ingredients. These excipients typically are inactive components that enhance the stability, solubility, and processing of drug formulations. Additionally, pharmaceutical excipients can impact the absorption, distribution, metabolism, and excretion (ADME) processes when drugs are administered together.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481