

Calcitriol-13C3

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Calcitriol-13C3 is a 13C-labeled form of calcitriol. Calcitriol (T6316) is the most potent metabolite of vitamin D and acts as an agonist of the vitamin D receptor (VDR).
Targets(IC50)	Sodium Channel

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481