

Tempol-D17,15N

Chemical Properties

CAS No. : 90429-66-6

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Tempol-D17,15N is a superoxide dismutase (SOD) mimic that effectively neutralizes reactive oxygen species (ROS).
Targets(IC50)	Reactive Oxygen Species, Autophagy

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481