


## Biotin-Ahx-Angiotensin II human

## Chemical Properties

CAS No. :	115102-72-2
Formula:	C66H96N16O15S
Molecular Weight:	1385.63
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Biotin-Ahx-Angiotensin II human (Biotin-Ahx-Angiotensin II; Biotin-Ahx-Ang II; Biotin-Ahx-DRVYIHPF) is a bioactive peptide.
Targets(IC50)	RAAS

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7217 mL	3.6085 mL	7.2169 mL
5 mM	0.1443 mL	0.7217 mL	1.4434 mL
10 mM	0.0722 mL	0.3608 mL	0.7217 mL
50 mM	0.0144 mL	0.0722 mL	0.1443 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481