

Tucidinostat-D4

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Tucidinostat-D4 is the deuterated form of Tucidinostat. Tucidinostat (T4481) is an effective, orally bioavailable inhibitor of HDAC class I HDAC1/2/3 and class IIb HDAC10, with IC50 values of 95, 160, 67, and 78 nM, respectively.
Targets(IC50)	HDAC

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481