

Caspofungin-D

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Caspofungin-D is the deuterated form of Caspofungin (T70350) diacetate. Caspofungin (T70350) (MK-0991) diacetate is a potent antifungal agent that inhibits the synthesis of 1,3-β-D-glucan synthase.
Targets(IC50)	Antibacterial, Antibiotic, Antifungal

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481