

(Rac)-Norepinephrine-D3 formate

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(Rac)-Norepinephrine-D3 formate is the deuterated form of Norepinephrine. Norepinephrine (Levarterenol; L-Noradrenaline) functions as a potent agonist for adrenergic receptors (AR). It can activate α 1, α 2, and β 1 receptors.
Targets(IC50)	Endogenous Metabolite,Adrenergic Receptor,Autophagy

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481