


## Devimistat-D10

## Chemical Properties

CAS No. :	2586055-61-8
Formula:	C <sub>22</sub> H <sub>28</sub> O <sub>2</sub> S <sub>2</sub>
Molecular Weight:	398.65
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Devimistat-D10 is the deuterium-labeled form of Devimistat (T6157).
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5085 mL	12.5423 mL	25.0847 mL
5 mM	0.5017 mL	2.5085 mL	5.0169 mL
10 mM	0.2508 mL	1.2542 mL	2.5085 mL
50 mM	0.0502 mL	0.2508 mL	0.5017 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481