

## Eribulin-D3 mesylate

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Eribulin-D3 mesylate is the deuterated form of Eribulin mesylate (T13687-GMP). Eribulin mesylate, an anti-cancer agent that targets microtubules (microtubule), is utilized for cancer research.
Targets(IC50)	Apoptosis, Microtubule Associated

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481