

## Idasanutlin (enantiomer)

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Idasanutlin enantiomer is an isomer of Idasanutlin and can be used as a control compound in experiments. Idasanutlin (RG7388) is a potent and selective MDM2 antagonist that inhibits the interaction between p53 and MDM2, with an IC50 value of 6 nM.
Targets(IC50)	Apoptosis,MDM-2/p53

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481