

DPPC-D9

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	DPPC-D9 is the deuterated form of DPPC. DPPC (T7773) (129Y83) is a phosphoglyceride used in the preparation of liposome monolayers. DPPC-liposomes serve effectively as delivery vectors for inducing immune responses against GSL antigens in mice.
Targets(IC50)	Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481