

## Danvatirsén sodium

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Danvantisen sodium is an antisense oligonucleotide targeting STAT3 with potential antineoplastic properties. It binds to STAT3 mRNA, thereby inhibiting the translation of the transcript. Suppressing STAT3 expression can lead to the induction of tumor cell apoptosis and inhibition of tumor cell growth.
Targets(IC50)	Apoptosis,STAT

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481