

(R,R)-S63845

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(R,R)-S63845 is the isomer of S63845 and serves as a control compound in experiments. S63845 is a potent and selective inhibitor of myeloid cell leukemia 1 (MCL1), binding to human MCL1 with a dissociation constant (Kd) of 0.19 nM.
Targets(IC50)	Bcl-2 Family

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481