


## HER2-IN-5

## Chemical Properties

CAS No. :	2691891-36-6
Formula:	C <sub>27</sub> H <sub>33</sub> N <sub>7</sub> O <sub>3</sub>
Molecular Weight:	503.6
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	HER2-IN-5 is an effective inhibitor of orally active HER-2.
Targets(IC50)	EGFR,Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9857 mL	9.9285 mL	19.857 mL
5 mM	0.3971 mL	1.9857 mL	3.9714 mL
10 mM	0.1986 mL	0.9929 mL	1.9857 mL
50 mM	0.0397 mL	0.1986 mL	0.3971 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481