

## V5 Epitope Tag Peptide Trifluoroacetate

## Chemical Properties

CAS No. :

Formula: C<sub>64</sub>H<sub>108</sub>N<sub>16</sub>O<sub>20</sub>.C<sub>2</sub>HF<sub>3</sub>O<sub>2</sub>

Molecular Weight: 1535.66

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

V5 Epitope Tag Peptide Trifluoroacetate is a Tag Peptide derived from a small Epitope on P and V proteins of monkey-virus paramyxovirus.

## Solubility Information

Solubility

H<sub>2</sub>O: Soluble,  
(< 1 mg/ml refers to the product slightly soluble or insoluble)

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6512 mL	3.2559 mL	6.5119 mL
5 mM	0.1302 mL	0.6512 mL	1.3024 mL
10 mM	0.0651 mL	0.3256 mL	0.6512 mL
50 mM	0.013 mL	0.0651 mL	0.1302 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481