

Ro 23-3423

## Chemical Properties


CAS No. : 110996-51-5

Formula: C<sub>23</sub>H<sub>22</sub>N<sub>4</sub>O<sub>2</sub>

Molecular Weight: 386.45

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ro 23-3423 is a specific inhibitor of thromboxane synthetase.
Targets(IC50)	Others, Prostaglandin Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5877 mL	12.9383 mL	25.8766 mL
5 mM	0.5175 mL	2.5877 mL	5.1753 mL
10 mM	0.2588 mL	1.2938 mL	2.5877 mL
50 mM	0.0518 mL	0.2588 mL	0.5175 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Morel DR, Huttemeier PC, Skoskiewicz MJ, Nguyenduy T, Melvin C, Robinson DR, Zapol WM. Dose-dependent effects of a pyridoquinazoline thromboxane synthetase inhibitor on arachidonic acid metabolites and hemodynamics during E. coli endotoxemia in anesthetized sheep. Prostaglandins. 1987 Jun;33(6):879-902. PubMed PMID: 3118413.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481